

## Anti-ERK 8 antibody



**Description** Rabbit polyclonal to ERK 8.

Model STJ92996

**Host** Rabbit

**Reactivity** Human

**Applications** ELISA, IF, IHC

**Immunogen** Synthesized peptide derived from human ERK 8

Immunogen Region 330-410 aa, Internal

**Gene ID** 225689

Gene Symbol MAPK15

**Dilution range** IHC 1:100-1:300IF 1:200-1:1000ELISA 1:5000

**Specificity** ERK 8 Polyclonal Antibody detects endogenous levels of ERK 8 protein.

**Tissue Specificity** Widely expressed with a maximal expression in lung and kidney.

**Purification** The antibody was affinity-purified from rabbit antiserum by affinity-

chromatography using epitope-specific immunogen.

**Note** For Research Use Only (RUO).

**Protein Name** Mitogen-activated protein kinase 15 MAP kinase 15 MAPK 15 Extracellular

signal-regulated kinase 7 ERK-7 Extracellular signal-regulated kinase 8

ERK-8

Molecular Weight 60 kDa

**Clonality** Polyclonal

**Conjugation** Unconjugated

**Isotype** IgG

**Formulation** Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

**Concentration** 1 mg/ml

**Storage Instruction** Store at -20°C, and avoid repeat freeze-thaw cycles.

Database Links <u>HGNC:24667</u>OMIM:NA

Alternative Names Mitogen-activated protein kinase 15 MAP kinase 15 MAPK 15 Extracellular

signal-regulated kinase 7 ERK-7 Extracellular signal-regulated kinase 8

ERK-8

**Function** In vitro, phosphorylates MBP.

**Sequence and Domain Family** The N-terminal region (1-20) is the minimal region necessary for

ubiquitination and further proteasomal degradation. The TXY motif contains the threonine and tyrosine residues whose phosphorylation activates the MAP

kinases.

**Post-translational** Dually phosphorylated on Thr-175 and Tyr-177, which activates the enzyme.

Autophosphorylated on threonine and tyrosine residues in vitro.

Ubiquitinated. Ubiquitination may allow its tight kinase activity regulation

and rapid turnover. May be ubiquitinated by a SCF E3 ligase.

St John's Laboratory Ltd

**Modifications** 

**F** +44 (0)207 681 2580

T +44 (0)208 223 3081

W http://www.stjohnslabs.com/ E info@stjohnslabs.com