


Anti-MMP-9 antibody





Description Rabbit polyclonal to MMP-9.

Model STJ93257

Host Rabbit

Reactivity Human

Applications ELISA, IF, IHC, WB

Immunogen Synthesized peptide derived from human MMP-9

Immunogen Region 620-700 aa, C-terminal

Gene ID 4318

Gene Symbol MMP9

Dilution range WB 1:500-1:2000IHC 1:100-1:300IF 1:200-1:1000ELISA 1:20000

Specificity MMP-9 Polyclonal Antibody detects endogenous levels of MMP-9 protein.

Tissue Specificity Produced by normal alveolar macrophages and granulocytes.

Purification The antibody was affinity-purified from rabbit antiserum by affinity-

chromatography using epitope-specific immunogen.

Note For Research Use Only (RUO).

Protein Name Matrix metalloproteinase-9 MMP-9 92 kDa gelatinase 92 kDa type IV

collagenase Gelatinase B GELB 67 kDa matrix metalloproteinase-9 82 kDa

matrix metalloproteinase-9

Molecular Weight 78 kDa

Clonality Polyclonal

Conjugation Unconjugated

Isotype IgG

Formulation Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Concentration 1 mg/ml

Storage Instruction Store at -20°C, and avoid repeat freeze-thaw cycles.

Database Links <u>HGNC:7176OMIM:120361</u>

Alternative Names Matrix metalloproteinase-9 MMP-9 92 kDa gelatinase 92 kDa type IV

collagenase Gelatinase B GELB 67 kDa matrix metalloproteinase-9 82 kDa

matrix metalloproteinase-9

Function May play an essential role in local proteolysis of the extracellular matrix and

in leukocyte migration. Could play a role in bone osteoclastic resorption. Cleaves KiSS1 at a Gly-|-Leu bond. Cleaves type IV and type V collagen into large C-terminal three quarter fragments and shorter N-terminal one quarter

fragments. Degrades fibronectin but not laminin or Pz-peptide.

Sequence and Domain Family The conserved cysteine present in the cysteine-switch motif binds the catalytic

zinc ion, thus inhibiting the enzyme. The dissociation of the cysteine from the

zinc ion upon the activation-peptide release activates the enzyme.

Cellular Localization Secreted, extracellular space, extracellular matrix

Post-translational Processing of the precursor yields different active forms of 64, 67 and 82 kDa.

Sequentially processing by MMP3 yields the 82 kDa matrix

metalloproteinase-9. N- and O-glycosylated.

St John's Laboratory Ltd

Modifications

F +44 (0)207 681 2580

T+44 (0)208 223 3081

W http://www.stjohnslabs.com/ E info@stjohnslabs.com