

## **Anti-HCN2 antibody**



**Description** Rabbit polyclonal to HCN2.

Model STJ93476

**Host** Rabbit

**Reactivity** Human, Mouse, Rat

**Applications** ELISA, IHC, WB

**Immunogen** Synthesized peptide derived from human HCN2

**Immunogen Region** 460-540 aa, Internal

**Gene ID** <u>610</u>

Gene Symbol HCN2

**Dilution range** WB 1:500-1:2000IHC 1:100-1:300ELISA 1:40000

**Specificity** HCN2 Polyclonal Antibody detects endogenous levels of HCN2 protein.

**Tissue Specificity** Highly expressed throughout the brain. Detected at low levels in heart.

**Purification** The antibody was affinity-purified from rabbit antiserum by affinity-

chromatography using epitope-specific immunogen.

**Note** For Research Use Only (RUO).

Protein Name Potassium/sodium hyperpolarization-activated cyclic nucleotide-gated channel

2 Brain cyclic nucleotide-gated channel 2 BCNG-2

Molecular Weight 100 kDa

**Clonality** Polyclonal

Conjugation Unconjugated

**Isotype IgG** 

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. **Formulation** 

1 mg/ml Concentration

Store at -20°C, and avoid repeat freeze-thaw cycles. **Storage Instruction** 

**Database Links** HGNC:4846OMIM:602781

Potassium/sodium hyperpolarization-activated cyclic nucleotide-gated channel **Alternative Names** 

2 Brain cyclic nucleotide-gated channel 2 BCNG-2

**Function** Hyperpolarization-activated ion channel exhibiting weak selectivity for

> potassium over sodium ions. Contributes to the native pacemaker currents in heart (If) and in neurons (Ih). Can also transport ammonium in the distal nephron. Produces a large instantaneous current. Modulated by intracellular chloride ions and pH; acidic pH shifts the activation to more negative voltages

**Sequence and Domain Family** The segment S4 is probably the voltage-sensor and is characterized by a series

of positively charged amino acids at every third position.

Cell membrane **Cellular Localization** 

Post-translational

Phosphorylation at Ser-668 by PRKG2 shifts the voltage-dependence to more negative voltages, hence counteracting the stimulatory effect of cGMP on **Modifications** 

gating.

St John's Laboratory Ltd

F +44 (0)207 681 2580

**T** +44 (0)208 223 3081

W http://www.stjohnslabs.com/ E info@stjohnslabs.com