

## Anti-Myosin XVI antibody



**Description** Rabbit polyclonal to Myosin XVI.

Model STJ94316

**Host** Rabbit

**Reactivity** Human

**Applications** ELISA, WB

Immunogen Synthesized peptide derived from human Myosin XVI

Immunogen Region 1050-1130 aa, Internal

**Gene ID** 23026

Gene Symbol MYO16

**Dilution range** WB 1:500-1:2000ELISA 1:10000

Specificity Myosin XVI Polyclonal Antibody detects endogenous levels of Myosin XVI

protein.

**Purification** The antibody was affinity-purified from rabbit antiserum by affinity-

chromatography using epitope-specific immunogen.

**Note** For Research Use Only (RUO).

Protein Name Unconventional myosin-XVI Neuronal tyrosine-phosphorylated

phosphoinositide-3-kinase adapter 3 Unconventional myosin-16

Molecular Weight 200 kDa

**Clonality** Polyclonal

**Conjugation** Unconjugated

**Isotype** IgG

**Formulation** Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

**Concentration** 1 mg/ml

**Storage Instruction** Store at -20°C, and avoid repeat freeze-thaw cycles.

Database Links <u>HGNC:29822OMIM:615479</u>

Alternative Names Unconventional myosin-XVI Neuronal tyrosine-phosphorylated

phosphoinositide-3-kinase adapter 3 Unconventional myosin-16

**Function** Myosins are actin-based motor molecules with ATPase activity.

Unconventional myosins serve in intracellular movements. Their highly divergent tails are presumed to bind to membranous compartments, which would be moved relative to actin filaments. May be involved in targeting of the catalytic subunit of protein phosphatase 1 during brain development. Activates PI3K and concomitantly recruits the WAVE1 complex to the close

vicinity of PI3K and regulates neuronal morphogenesis.

**Cellular Localization** Cytoplasm. Found in puncta in soma and processes of astrocytes and

dissociated cerebellar cells with the morphology of migrating granule cells.

**Post-translational** 

Modifications

Phosphorylated on tyrosine residues by FYN upon stimulation with CNTN5.

St John's Laboratory Ltd

**F** +44 (0)207 681 2580

T+44 (0)208 223 3081

W http://www.stjohnslabs.com/ E info@stjohnslabs.com