


## Anti-Rabenosyn-5 antibody



**Description** Rabbit polyclonal to Rabenosyn-5.

Model STJ95314

**Host** Rabbit

**Reactivity** Human, Mouse

**Applications** ELISA, WB

Immunogen Synthesized peptide derived from human Rabenosyn-5

Immunogen Region 120-200 aa, Internal

**Gene ID** <u>64145</u>

Gene Symbol RBSN

**Dilution range** WB 1:500-1:2000ELISA 1:10000

**Specificity** Rabenosyn-5 Polyclonal Antibody detects endogenous levels of Rabenosyn-5

protein.

**Purification** The antibody was affinity-purified from rabbit antiserum by affinity-

chromatography using epitope-specific immunogen.

**Note** For Research Use Only (RUO).

**Protein Name** Rabenosyn-5 110 kDa protein FYVE finger-containing Rab5 effector protein

rabenosyn-5 RAB effector RBSN Zinc finger FYVE domain-containing

protein 20

Molecular Weight 89 kDa

**Clonality** Polyclonal

**Conjugation** Unconjugated

**Isotype** IgG

**Formulation** Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

**Concentration** 1 mg/ml

**Storage Instruction** Store at -20°C, and avoid repeat freeze-thaw cycles.

Database Links <u>HGNC:20759OMIM:609511</u>

Alternative Names Rabenosyn-5 110 kDa protein FYVE finger-containing Rab5 effector protein

rabenosyn-5 RAB effector RBSN Zinc finger FYVE domain-containing

protein 20

**Function** Rab4/Rab5 effector protein acting in early endocytic membrane fusion and

membrane trafficking of recycling endosomes. Required for endosome fusion either homotypically or with clathrin coated vesicles. Plays a role in the

lysosomal trafficking of CTSD/cathepsin D from the Golgi to lysosomes. Also promotes the recycling of transferrin directly from early endosomes to the

plasma membrane. Binds phospholipid vesicles containing

phosphatidylinositol 3-phosphate (PtdInsP3) . Plays a role in the recycling of

transferrin receptor to the plasma membrane.

Cellular Localization Cell membrane. Lipid-anchor. Cytoplasmic side. Early endosome membrane

St John's Laboratory Ltd

**F** +44 (0)207 681 2580

T+44 (0)208 223 3081

W http://www.stjohnslabs.com/ E info@stjohnslabs.com