

## **Anti-SENP6** antibody



**Description** Rabbit polyclonal to SENP6.

Model STJ95605

**Host** Rabbit

**Reactivity** Human

**Applications** ELISA, IHC

ImmunogenSynthesized peptide derived from human SENP6

Immunogen Region 1010-1090 aa, C-terminal

**Gene ID** <u>26054</u>

Gene Symbol SENP6

**Dilution range** IHC 1:100-1:300ELISA 1:10000

**Specificity** SENP6 Polyclonal Antibody detects endogenous levels of SENP6 protein.

**Tissue Specificity** Highly expressed in reproductive organs, such as testis, ovary and prostate.

**Purification** The antibody was affinity-purified from rabbit antiserum by affinity-

chromatography using epitope-specific immunogen.

**Note** For Research Use Only (RUO).

Protein Name Sentrin-specific protease 6 SUMO-1-specific protease 1 Sentrin/SUMO-

specific protease SENP6

Molecular Weight 126.146 kDa

**Clonality** Polyclonal

**Conjugation** Unconjugated

**Isotype** IgG

**Formulation** Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

**Concentration** 1 mg/ml

**Storage Instruction** Store at -20°C, and avoid repeat freeze-thaw cycles.

Database Links <u>HGNC:20944OMIM:605003</u>

Alternative Names Sentrin-specific protease 6 SUMO-1-specific protease 1 Sentrin/SUMO-

specific protease SENP6

**Function** Protease that deconjugates SUMO1, SUMO2 and SUMO3 from targeted

proteins. Processes preferentially poly-SUMO2 and poly-SUMO3 chains, but does not efficiently process SUMO1, SUMO2 and SUMO3 precursors. Deconjugates SUMO1 from RXRA, leading to transcriptional activation. Involved in chromosome alignment and spindle assembly, by regulating the kinetochore CENPH-CENPI-CENPK complex. Desumoylates PML and CENPI, protecting them from degradation by the ubiquitin ligase RNF4, which targets polysumoylated proteins for proteasomal degradation.

Desumoylates also RPA1, thus preventing recruitment of RAD51 to the DNA  $\,$ 

damage foci to initiate DNA repair through homologous recombination.

**Cellular Localization** Nucleus

St John's Laboratory Ltd

F +44 (0)207 681 2580

**T** +44 (0)208 223 3081

W http://www.stjohnslabs.com/ E info@stjohnslabs.com