


Anti-UBA2 antibody





Description Rabbit polyclonal to UBA2.

Model STJ96159

Host Rabbit

Reactivity Human

Applications ELISA, IF, IHC, WB

Immunogen Synthesized peptide derived from human UBA2

Immunogen Region 560-640 aa, C-terminal

Gene ID <u>10054</u>

Gene Symbol <u>UBA2</u>

Dilution range WB 1:500-1:2000IHC 1:100-1:300IF 1:200-1:1000ELISA 1:10000

Specificity UBA2 Polyclonal Antibody detects endogenous levels of UBA2 protein.

Purification The antibody was affinity-purified from rabbit antiserum by affinity-

chromatography using epitope-specific immunogen.

Note For Research Use Only (RUO).

Protein Name SUMO-activating enzyme subunit 2 Anthracycline-associated resistance ARX

Ubiquitin-like 1-activating enzyme E1B Ubiquitin-like modifier-activating

enzyme 2

Molecular Weight 71 kDa

Clonality Polyclonal

Conjugation Unconjugated

Isotype IgG

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. **Formulation**

1 mg/ml Concentration

Store at -20°C, and avoid repeat freeze-thaw cycles. **Storage Instruction**

Database Links HGNC:30661OMIM:613295

Alternative Names SUMO-activating enzyme subunit 2 Anthracycline-associated resistance ARX

Ubiquitin-like 1-activating enzyme E1B Ubiquitin-like modifier-activating

enzyme 2

The heterodimer acts as an E1 ligase for SUMO1, SUMO2, SUMO3, and **Function**

> probably SUMO4. It mediates ATP-dependent activation of SUMO proteins followed by formation of a thioester bond between a SUMO protein and a

conserved active site cysteine residue on UBA2/SAE2.

Cellular Localization Cytoplasm. Nucleus. Shuttles between the cytoplasm and the nucleus,

sumoylation is required either for nuclear translocation or nuclear retention.

Sumoylated with SUMO1 and SUMO2/3 and by UBC9. Sumoylation at Post-translational

Lys-236 inhibits enzymatic activity. Sumoylation at the C-terminal lysine **Modifications**

cluster plays an essential role in nuclear trafficking.

St John's Laboratory Ltd

F +44 (0)207 681 2580

W http://www.stjohnslabs.com/ T+44 (0)208 223 3081 E info@stjohnslabs.com