

Anti-TAF II p140 antibody



Description Rabbit polyclonal to TAF II p140.

Model STJ96387

Host Rabbit

Human, Mouse, Rat Reactivity

ELISA, WB **Applications**

Synthesized peptide derived from human TAF II p140. **Immunogen**

Internal **Immunogen Region**

Gene ID 83860

Gene Symbol TAF3

Dilution range WB 1:500-1:2000ELISA 1:10000

TAF II p140 Polyclonal Antibody detects endogenous levels of TAF II p140 **Specificity**

protein.

The antibody was affinity-purified from rabbit antiserum by affinity-**Purification**

chromatography using epitope-specific immunogen.

For Research Use Only (RUO). Note

Protein Name Transcription initiation factor TFIID subunit 3 140 kDa TATA box-binding

> protein-associated factor TBP-associated factor 3 Transcription initiation factor TFIID 140 kDa subunit TAF II140 TAF140 TAFII-140 TAF

100 kDa **Molecular Weight**

Polyclonal **Clonality**

Conjugation Unconjugated

Isotype IgG

Formulation Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Concentration 1 mg/ml

Storage Instruction Store at -20°C, and avoid repeat freeze-thaw cycles.

Database Links <u>HGNC:17303OMIM:606576</u>

Alternative Names Transcription initiation factor TFIID subunit 3 140 kDa TATA box-binding

protein-associated factor TBP-associated factor 3 Transcription initiation factor TFIID 140 kDa subunit TAF II140 TAF140 TAFII-140 TAF

Function Transcription factor TFIID is one of the general factors required for accurate

and regulated initiation by RNA polymerase II. TFIID is a multimeric protein complex that plays a central role in mediating promoter responses to various

activators and repressors. Required in complex with TBPL2 for the

differentiation of myoblasts into myocytes. The complex replaces TFIID at

specific promoters at an early stage in the differentiation process.

Sequence and Domain Family The PHD-type zinc finger mediates binding to histone H3 methyllysine at

position 4 (H3K4me3).

Cellular Localization Nucleus

St John's Laboratory Ltd

F +44 (0)207 681 2580 **T** +44 (0)208 223 3081

W http://www.stjohnslabs.com/ E info@stjohnslabs.com