


## Anti-CyPA antibody





Description	Rabbit polyclonal to CyPA.

Model STJ96562

**Host** Rabbit

**Reactivity** Human, Mouse, Rat

**Applications** ELISA, IHC, WB

**Immunogen** Synthesized peptide derived from human CyPA.

**Immunogen Region** 51-100 aa, Internal

**Gene ID** <u>5478</u>

Gene Symbol PPIA

**Dilution range** WB 1:500-1:2000IHC-P 1:100-300ELISA 1:20000

**Specificity** CyPA Polyclonal Antibody detects endogenous levels of CyPA protein.

**Purification** The antibody was affinity-purified from rabbit antiserum by affinity-

chromatography using epitope-specific immunogen.

**Note** For Research Use Only (RUO).

Protein Name Peptidyl-prolyl cis-trans isomerase A PPIase A Cyclophilin A Cyclosporin A-

binding protein Rotamase A Peptidyl-prolyl cis-trans isomerase A, N-

terminally processed

Molecular Weight 18 kDa

**Clonality** Polyclonal

**Conjugation** Unconjugated

**Isotype** IgG

**Formulation** Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

**Concentration** 1 mg/ml

**Storage Instruction** Store at -20°C, and avoid repeat freeze-thaw cycles.

Database Links <u>HGNC:9253OMIM:123840</u>

Alternative Names Peptidyl-prolyl cis-trans isomerase A PPIase A Cyclophilin A Cyclosporin A-

binding protein Rotamase A Peptidyl-prolyl cis-trans isomerase A, N-

terminally processed

**Function** PPIases accelerate the folding of proteins. It catalyzes the cis-trans

isomerization of proline imidic peptide bonds in oligopeptides.

**Cellular Localization** Cytoplasm Secreted. Secretion occurs in response to oxidative stress in

vascular smooth muscle through a vesicular secretory pathway that involves actin remodeling and myosin II activation, and mediates ERK1/2 activation.

**Post-translational** Acetylation at Lys-125 markedly inhibits catalysis of cis to trans

**Modifications** isomerization and stabilizes cis rather than trans forms of the HIV-1 capsid.

PPIA acetylation also antagonizes the immunosuppressive effects of cyclosporine by inhibiting the sequential steps of cyclosporine binding and

calcineurin inhibition.

St John's Laboratory Ltd

**F** +44 (0)207 681 2580

T+44 (0)208 223 3081

W http://www.stjohnslabs.com/ E info@stjohnslabs.com