


## Anti-MMP-9 antibody





**Description** Mouse monoclonal to MMP-9.

Model STJ98249

**Host** Mouse

**Reactivity** Human

**Applications** ELISA, FC, WB

**Immunogen** Purified recombinant fragment of human MMP-9 expressed in E. Coli.

**Gene ID** 4318

Gene Symbol MMP9

**Dilution range** WB 1:500-1:2000FC 1:200-1:400ELISA 1:10000

**Specificity** MMP-9 Monoclonal Antibody detects endogenous levels of MMP-9 protein.

**Tissue Specificity** Produced by normal alveolar macrophages and granulocytes.

**Purification** Affinity purification

Clone ID 5C3

**Note** For Research Use Only (RUO).

**Protein Name** Matrix metalloproteinase-9 MMP-9 92 kDa gelatinase 92 kDa type IV

collagenase Gelatinase B GELB 67 kDa matrix metalloproteinase-9 82 kDa

matrix metalloproteinase-9

**Clonality** Monoclonal

**Conjugation** Unconjugated

**Isotype** IgG2a

**Formulation** Ascitic fluid containing 0.03% sodium azide.

**Storage Instruction** Store at -20°C, and avoid repeat freeze-thaw cycles.

Database Links <u>HGNC:71760MIM:120361</u>

Alternative Names Matrix metalloproteinase-9 MMP-9 92 kDa gelatinase 92 kDa type IV

collagenase Gelatinase B GELB 67 kDa matrix metalloproteinase-9 82 kDa

matrix metalloproteinase-9

**Function** May play an essential role in local proteolysis of the extracellular matrix and

in leukocyte migration. Could play a role in bone osteoclastic resorption. Cleaves KiSS1 at a Gly-|-Leu bond. Cleaves type IV and type V collagen into large C-terminal three quarter fragments and shorter N-terminal one quarter

fragments. Degrades fibronectin but not laminin or Pz-peptide.

Sequence and Domain Family The conserved cysteine present in the cysteine-switch motif binds the catalytic

zinc ion, thus inhibiting the enzyme. The dissociation of the cysteine from the

zinc ion upon the activation-peptide release activates the enzyme.

Cellular Localization Secreted, extracellular space, extracellular matrix

**Post-translational** Processing of the precursor yields different active forms of 64, 67 and 82 kDa.

**Modifications** Sequentially processing by MMP3 yields the 82 kDa matrix

metalloproteinase-9. N- and O-glycosylated.

St John's Laboratory Ltd

F+44 (0)207 681 2580

**T** +44 (0)208 223 3081

W http://www.stjohnslabs.com/ E info@stjohnslabs.com