# Data Sheet (Cat.No.T1843)



## XMD8-92

## **Chemical Properties**


CAS No.: 1234480-50-2

Formula: C26H30N6O3

Molecular Weight: 474.55

Appearance: no data available

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



## **Biological Description**

Description	XMD8-92 is an effective and specific BMK1/ERK5 inhibitor (Kd: 80 nM).
Targets(IC50)	ERK,BMI-1,Epigenetic Reader Domain
In vitro	XMD8-92 significantly inhibits the growth of pancreatic tumor xenografts by downregulating DCLK1 and several of its downstream targets. At a dose of 50 mg/kg, administered intraperitoneally (i.p.), XMD8-92 blocks tumor cell proliferation and tumorassociated angiogenesis, resulting in significant suppression of tumor growth in both xenografted human and syngeneic mouse models.
In vivo	Hydroxysafflor yellow A (HSYA) inhibits the activation of hepatic stellate cells, while XMD8-92 can significantly block this inhibition and prevent the HSYA-mediated downregulation of MEF2C. By inhibiting the activation of BMK1, XMD8-92 significantly induces the expression of p21 in cells and inhibits the proliferation of cancer cells.
Kinase Assay	KiNativ profiling of XMD8-92 is carried out with both an ATP and ADP acylphosphate-desthiobiotin with the following modifications. HeLa cell lysates (5 mg/mL total protein) are incubated in the presence of XMD8-92 at 50 μM, 10 μM, 2 μM, 0.8 μM, and 0 μM for 15 minutes prior to addition of the ATP or ADP acylphosphate probe (5 μM final probe concentration). All reactions are performed in duplicate. Probe reactions proceeded for 10 minutes and the reaction stopped by the addition of urea and processed for MS analysis. Samples are analyzed by LC-MS/MS on a linear ion trap mass spectrometer using a time segmented "target list" designed to collect MS/MS spectra from all kinase peptide-probe conjugates that can be detected in HeLa cell lysates. This target list is generated and validated by prior exhaustive analysis of HeLa lysates. Up to four characteristic fragment ions for each kinase peptide-probe conjugate are used to extract signals for each kinase, and a comparison of inhibitor treated to control (untreated) lysates allow for precise determination of % inhibition at each point. A manuscript describing the details of this targeted mass spectrometry approach is in preparation[1].

## **Solubility Information**

#### A DRUG SCREENING EXPERT

	Solubility	DMSO: 19 mg/mL (40.04 mM), Sonication is recommended.
		(< 1 mg/ml refers to the product slightly soluble or insoluble)

#### **Preparing Stock Solutions**

	1mg	5mg	10mg	
1 mM	2.1073 mL	10.5363 mL	21.0726 mL	
5 mM	0.4215 mL	2.1073 mL	4.2145 mL	
10 mM	0.2107 mL	1.0536 mL	2.1073 mL	
50 mM	0.0421 mL	0.2107 mL	0.4215 mL	

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

## Reference

Yang Q, et al. Cancer Cell. 2010, 18(3), 258-267. Dong H, et al. Pharm Biol. 2013. doi:10.3109/13880209. Sureban SM, et al. Cancer Lett. 2014, 351(1), 151-161.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481

Page 2 of 2 www.targetmol.com