


## SB-3CT

## Chemical Properties

CAS No. : 292605-14-2  
 Formula: C<sub>15</sub>H<sub>14</sub>O<sub>3</sub>S<sub>2</sub>  
 Molecular Weight: 306.4  
 Appearance: no data available  
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



## Biological Description

Description	SB-3CT is a potent and specific gelatinase inhibitor for MMP-2 (K <sub>i</sub> = 13.9 nM) and MMP-9 (K <sub>i</sub> = 600 nM).
Targets(IC <sub>50</sub> )	MMP
In vitro	In a permanent cerebral ischemia mouse model induced by an embolus, SB-3CT protects neurons from ischemic cell death by inhibiting the degradation of neuronal adhesion proteins. Furthermore, daily intraperitoneal injections of SB-3CT (5-50 mg/kg) in the L-CI.5sT murine lymphoma model significantly hinders liver cancer cell metastasis. Likewise, in a model of osteolytic metastasis of human prostate cancer cells, daily intraperitoneal administration of SB-3CT (50 mg/kg) curtails cell growth, bone resorption, and angiogenesis.
In vivo	SB-3CT is capable of reducing osteolytic responses, aiding in the maintenance of bone integrity. In PC3 tumor cells, SB-3CT inhibits cellular growth. Within the basement membrane matrix, SB-3CT directly suppresses the infiltration of bone marrow endothelial cells and tubule formation.
Kinase Assay	kinase assays: PI3K inhibition by LY294002 is determined in a radiometric assay using purified, recombinant enzymes with 1 μM ATP. The kinase reaction is carried out for 1 hour at room temperature (24°C) and is terminated by addition of PBS. IC <sub>50</sub> values are subsequently determined using a sigmoidal dose-response curve fit (variable slope). CK2 and GSK3β (glycogen synthase kinase 3β) inhibition is established by kinase selectivity screening. LY294002 is tested against the Upstate panel of kinases in 10 μM ATP.
Cell Research	SB-3CT is dissolved in DMSO (10 mM) and stored (-22°C), and then diluted with appropriate media (DMSO 1%) before use[2]. PC3 cells are seeded in 35-mm dishes (5×10 <sup>4</sup> cells/dish) in complete culture medium. The next day, the medium is replaced with complete medium supplemented with 1% DMSO alone (vehicle) or SB-3CT (final concentrations 0.1-50 μM) in 1% DMSO. At various times, the cells are harvested with trypsin and counted[2].

## Solubility Information

## A DRUG SCREENING EXPERT

Solubility	Ethanol: 6.1 mg/mL (19.91 mM), Sonication is recommended. DMSO: 55 mg/mL (179.5 mM), Sonication is recommended. (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	---

### Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.2637 mL	16.3185 mL	32.6371 mL
5 mM	0.6527 mL	3.2637 mL	6.5274 mL
10 mM	0.3264 mL	1.6319 mL	3.2637 mL
50 mM	0.0653 mL	0.3264 mL	0.6527 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

### Reference

Brown S, et al. J Am Chem Soc. 2000, 122 (28), 6799-6800.

Lei Z, Chen H, Zhang H, et al. Evaluation of Matrix Metalloproteinase Inhibition by Peptide Microarray-Based Fluorescence Assay on Polymer Brush Substrate and in Vivo Assessment. ACS applied materials & interfaces. 2017 Dec 20;9(50):44241-44250.

Lei Z, Jian M, Wei J, et al. Array-Based in situ Fluorescence Assay for Profiling Multiplex Matrix Metalloproteinases Activities in Tissue Section. Analytica Chimica Acta. 2019

Krüger A, et al. Cancer Res. 2005, 65(9), 3523-3526

Bonfil RD, et al. Int J Cancer. 2006, 118(11), 2721-2726.

Liu H, Yao C, Zhang L, et al. Nanoliposomes co-encapsulating Ce6 and SB3CT against the proliferation and metastasis of melanoma with the integration of photodynamic therapy and NKG2D-related immunotherapy on A375 cells. Nanotechnology. 2021

Cui J, et al. Mol Neurodegener. 2012, 15, 7-21.

Lei Z, Jian M, Wei J, et al. Array-Based in situ Fluorescence Assay for Profiling Multiplex Matrix Metalloproteinases Activities in Tissue Section[J]. Analytica Chimica Acta. 2019.

Lei Z, Chen H, Zhang H, et al. Evaluation of Matrix Metalloproteinase Inhibition by Peptide Microarray-Based Fluorescence Assay on Polymer Brush Substrate and in Vivo Assessment[J]. ACS applied materials & interfaces. 2017 Dec 20;9(50):44241-44250.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

**This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use**

Tel: 781-999-4286    E\_mail: info@targetmol.com    Address: 36 Washington Street, Wellesley Hills, MA 02481