


Talazoparib

Chemical Properties

CAS No. :	1207456-01-6
Formula:	C ₁₉ H ₁₄ F ₂ N ₆ O
Molecular Weight:	380.35
Appearance:	no data available
Storage:	keep away from moisture, store at low temperature Powder: -20°C for 3 years In solvent: -80°C for 1 year



Biological Description

Description	Talazoparib (LT-673) is a new-type PARP inhibitor (IC ₅₀ : 0.58 nM), It similarly binds to PARP1/2 (K _i : 1.2/0.85 nM).
Targets(IC ₅₀)	PARP
In vitro	Talazoparib (BMN 673) demonstrates excellent potency, inhibiting PARP1 and PARP2 enzyme activity with K _i = 1.2 and 0.87 nM, respectively. It inhibits PARP-mediated PARylation in a whole-cell assay with an EC ₅₀ of 2.51 nM and prevents the proliferation of cancer cells carrying mutant BRCA1/2, with EC ₅₀ = 0.3 nM (MX-1) and 5 nM (Capan-1), respectively [1]. BMN 673 is a potent PARP1/2 inhibitor (PARP1 IC ₅₀ = 0.57 nmol/L), but it does not inhibit other enzymes that we have tested. BMN 673 selectively targeted tumor cells with BRCA1, BRCA2, or PTEN gene defects with 20- to more than 200-fold greater potency than existing PARP1/2 inhibitors [2].
In vivo	Talazoparib is orally available, displaying favorable pharmacokinetic (PK) properties and remarkable antitumor efficacy in the BRCA1 mutant MX-1 breast cancer xenograft model following oral administration as a single agent or in combination with chemotherapy agents such as temozolomide and cisplatin [1]. Oral administration of BMN 673 elicited remarkable antitumor activity in vivo; xenografted tumors that carry defects in DNA repair due to BRCA mutations or PTEN deficiency were profoundly sensitive to oral BMN 673 treatment at well-tolerated doses in mice. Synergistic or additive antitumor effects were also found when BMN 673 was combined with temozolomide, SN38, or platinum drugs [2].
Kinase Assay	The ability of a test compound to inhibit PARP1 enzyme activity was assessed using the PARP Assay Kit following the manufacturer's instruction. IC ₅₀ values were calculated using GraphPad Prism5 software. For PARP inhibitor K _i determination, enzyme assays were conducted in 96-well FlashPlate with 0.5 U PARP1 enzyme, 0.25× activated DNA, 0.2 µCi [3H] NAD, and 5 µmol/L cold NAD in a final volume of 50 µL reaction buffer containing 10% glycerol (v/v), 25 mmol/L HEPES, 12.5 mmol/L MgCl ₂ , 50 mmol/L KCl, 1 mmol/L dithiothreitol (DTT), and 0.01% NP-40 (v/v), pH 7.6. Reactions were initiated by adding NAD to the PARP reaction mixture with or without inhibitors and incubated for 1 minute at room temperature. Fifty microliter of ice-cold 20% trichloroacetic acid (TCA) was then added to each well to stop the reaction. The plate was sealed and shaken for a further 120 minutes at room temperature, followed by centrifugation. Radioactive signal bound to the FlashPlate was determined using TopCount. PARP1 K _m was determined

using Michaelis-Menten equation from various substrate concentrations (1–100 $\mu\text{mol/L}$ NAD). Compound K_i was calculated from enzyme inhibition curve according to the formula: $K_i = \text{IC}_{50}/[1+(\text{substrate})/K_m]$. K_m for PARP2 enzyme and compound K_i were determined with the same assay protocol except 30 ng PARP2, 0.25 \times activated DNA, 0.2 μCi $[^3\text{H}]$ NAD, and 20 $\mu\text{mol/L}$ cold NAD were used in the reaction for 30 minutes at room temperature [2].

Cell Research	Colony formation assays were conducted as described previously. In brief, cells were seeded into 6-well plates at a concentration of 500 to 2,000 cells per well. After 24 hours, media was replaced with fresh media containing PARP1/2 inhibitor. This procedure was repeated twice weekly for 14 days, at which point colonies were fixed with TCA and stained with sulforhodamine B. Colonies were counted and surviving fractions calculated by normalizing colony counts to colony numbers in vehicle-treated wells. Survival curves were plotted using a four-parameter logistic regression curve fit [2].
Animal Research	Female athymic nu/nu mice (8–10-week old) were used for all in vivo xenograft studies. Mice were quarantined for at least 1 week before experimental manipulation. Exponentially growing cells (LNCap and MDA-MB-468) or in vivo passaged tumor fragments (MX-1) were implanted subcutaneously at the right flank of nude mice. When tumors reached an average volume of approximately 150 mm^3 , mice were randomized into various treatment groups (6–8 mice/group) in each study. Mice were visually observed daily and tumors were measured twice weekly by calliper to determine tumor volume using the formula $[\text{length}/2] \times [\text{width}^2]$. Group median tumor volume (mm^3) was graphed over time to monitor tumor growth. In single-agent studies, olaparib (100 mg/kg), BMN 673 (various doses as indicated), or vehicle (10% DMSO, 6% Solutol, and 84% PBS) was administered by oral gavage (per os), once daily or BMN 673 (0.165 mg/kg) twice daily for 28 consecutive days. Mice were continuously monitored for 10 more days after last day of dosing. In cisplatin combination study, BMN 673, olaparib, or vehicle was administered per os once daily for 8 days starting on day 1. Cisplatin at a dosage of 6 mg/kg or its vehicle (saline) was administered intraperitoneally as a single injection on day 3, 30 minutes after PARP inhibitor was administered. Combination with carboplatin was conducted in a similar way in MX-1 model in which BMN 673 was administered per os once daily for either 8 days or 5 days and carboplatin was injected intraperitoneally at single dose of 35 mg/kg, 30 minutes after BMN 673 on day 3 [2].

Solubility Information

Solubility	DMSO: 36 mg/mL (94.65 mM), Sonication is recommended. H ₂ O: < 1 mg/mL (insoluble or slightly soluble), Ethanol: < 1 mg/mL (insoluble or slightly soluble), 10% DMSO+90% Saline: 0.1 mg/mL (0.26 mM), Solution. (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.6292 mL	13.1458 mL	26.2916 mL
5 mM	0.5258 mL	2.6292 mL	5.2583 mL
10 mM	0.2629 mL	1.3146 mL	2.6292 mL
50 mM	0.0526 mL	0.2629 mL	0.5258 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Reference

- Wilkerson PM, et al. Preclinical evaluation of the PARP inhibitor BMN-673 for the treatment of ovarian clear cell cancer. *Oncotarget*. 2017 Jan 24;8(4):6057-6066.
- Hou X, Tian M, Ning J, et al. PARP inhibitor shuts down the global translation of thyroid cancer through promoting Pol II binding to DIMT1 pause. *International Journal of Biological Sciences*. 2023, 19(12): 3970.
- Shen Y, et al. BMN 673, a novel and highly potent PARP1/2 inhibitor for the treatment of human cancers with DNA repair deficiency. *Clin Cancer Res*. 2013 Sep 15;19(18):5003-15.
- Mateos-Pujante A, Jiménez M C, Andreu I. Assessment of the PARP inhibitor talazoparib photosafety profile. *Biomedicine & Pharmacotherapy*. 2023, 167: 115593.
- Kamaletdinova T, Zong W, Urbánek P, et al. Poly (ADP-Ribose) Polymerase-1 Lacking Enzymatic Activity Is Not Compatible with Mouse Development. *Cells*. 2023, 12(16): 2078.
- Kiehlmeier S, Rafiee M R, Bakr A, et al. Identification of therapeutic targets of the hijacked super-enhancer complex in EVI1-rearranged leukemia[J]. *Leukemia*. 2021: 1-12.
- Kiehlmeier S, Rafiee M R, Bakr A, et al. Identification of therapeutic targets of the hijacked super-enhancer complex in EVI1-rearranged leukemia. *Leukemia*. 2021 Nov;35(11):3127-3138. doi: 10.1038/s41375-021-01235-z. Epub 2021 Apr 28.
- Liu C, Li J, Xu F, et al. PARP1-DOT1L transcription axis drives acquired resistance to PARP inhibitor in ovarian cancer. *Molecular Cancer*. 2024, 23(1): 111.
- Liu S, Deng P, Yu Z, et al. CDC7 Inhibition Potentiates Antitumor Efficacy of PARP Inhibitor in Advanced Ovarian Cancer. *Advanced Science*. 2024: 2403782.
- Lodovichi S, Nepomuceno T C, Woods N T, et al. SART1 modulates Poly-(ADP-Ribose) chain accumulation and PARP1 chromatin localization. *iScience*. 2024

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 36 Washington Street, Wellesley Hills, MA 02481