


CCT128930

## Chemical Properties

CAS No. : 885499-61-6  
 Formula: C<sub>18</sub>H<sub>20</sub>ClN<sub>5</sub>  
 Molecular Weight: 341.84  
 Appearance: no data available  
 Storage: store at low temperature  
 Powder: -20°C for 3 years | In solvent: -80°C for 1 year



## Biological Description

Description	CCT128930 is a potent, ATP-competitive and selective inhibitor of Akt2 with IC <sub>50</sub> of 6 nM, 28-fold greater selectivity for Akt2 than the closely related PKA kinase.
Targets(IC <sub>50</sub> )	Akt, Autophagy, PKA, S6 Kinase
In vitro	CCT128930 exhibits marked antiproliferative activity against PTEN-deficient human tumor cell lines including U87 mg human glioblastoma cells, LNCaP human prostate cancer cells and PC3 human prostate cancer cells with GI <sub>50</sub> of 6.3 μM, 0.35 μM and 1.9 μM, respectively. Furthermore, CCT128930 causes a G1 arrest in PTEN-null U87 mg human glioblastoma cells and Akt pathway blockade. [1]
In vivo	CCT128930 at 25 mg/kg i.p. shows a marked antitumor effect in established PTEN-null U87 mg human glioblastoma xenografts with a treated:control (T/C) ratio of 48% on day 12. In HER2-positive, PIK3CA-mutant BT474 human breast cancer xenografts, CCT128930 at 40 mg/kg also produces a profound antitumor effect with complete growth arrest and a T/C ratio of 29% on day 22. CCT128930 administered via i.v. reaches a peak concentration of 6.4 μM in plasma and is eliminated with a relatively short half-life, high volume of distribution, and rapid clearance, giving an area under the curve AUC <sub>0-∞</sub> of 4.6 μM·h. CCT128930 administered via i.p. leads to the peak plasma drug concentration of 1.3 μM and the corresponding AUC <sub>0-∞</sub> of 1.3 μM·h. Oral CCT128930 administration leads to the peak plasma concentration of only 0.43 μM and a correspondingly low AUC <sub>0-∞</sub> of 0.4 μM·h. [1]
Kinase Assay	Kinase assays: Profiling against 50 different human kinases is carried out using 10 μM CCT128930 at an ATP concentration equivalent to the K <sub>m</sub> for each enzyme.
Cell Research	Cells are seeded in 96-well plates and allowed to attach for 36 hours to ensure exponential growth prior to treatment. In vitro antiproliferative activity is determined using a 96-hour SRB assay. TCA-fixed cells are stained for 30 minutes with 0.4% (wt/vol) SRB dissolved in 1% acetic acid. At the end of the staining period, SRB is removed and cultures are quickly rinsed four times with 1% acetic acid to remove unbound dye. The acetic acid is poured directly into the culture wells from a beaker. This procedure permits rinsing to be performed quickly so that desorption of protein-bound dye does not occur. Residual wash solution is removed by sharply flicking plates over a sink, which ensures the complete removal of rinsing solution. Because of the strong capillary action in 96-well plates, draining by gravity alone often fails to remove the rinse

solution when plates are simply inverted. After being rinsed, the cultures are air dried until no standing moisture is visible. Bound dye is solubilized with 10 mM unbuffered Tris base (pH 10.5) for 5 minutes on a gyratory shaker. OD is read in either a UVmax microtiter plate reader or a Beckman DU-70 spectrophotometer. For maximum sensitivity, OD is measured at 564 nm. Because readings are linear with dye concentrations only below 1.8 OD units, however, suboptimal wavelengths are generally used, so that all samples in an experiment remains within the linear OD range. With most cell lines, wavelengths of approximately 490-530 nm works well for this purpose.(Only for Reference)

### Solubility Information

Solubility	Ethanol: 6 mg/mL (17.55 mM),Sonication is recommended. H2O: <1 mg/mL, DMSO: 25 mg/mL (73.13 mM),Sonication is recommended. (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

### Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.9253 mL	14.6267 mL	29.2535 mL
5 mM	0.5851 mL	2.9253 mL	5.8507 mL
10 mM	0.2925 mL	1.4627 mL	2.9253 mL
50 mM	0.0585 mL	0.2925 mL	0.5851 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

### Reference

Yap TA, et al. Mol Cancer Ther. 2011, 10(2), 360-371.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481