# Data Sheet (Cat.No.TMPY-03180)



# Ribonuclease A Protein, Human, Recombinant (His)

#### **General Information**

Synonyms: RNS1;ribonuclease, RNase A family, 1 (pancreatic);RIB1;RAC1

Protein Construction:

A DNA sequence encoding the human RNASE1 (P07998) (Met1-Thr156) was expressed with a

polyhistidine tag at the C-terminus. Predicted N terminal: Lys 29

Species: Human

Expression Host: HEK293 Cells

Accession: P07998

Molecular Weight: 16 kDa (predicted); 29,26 and 22 kDa (reducing conditions)

## **QC** Testing

Measured by its ability to hydrolyze yeast RNA. One unit of the enzyme causes an increase in

Biological Activity: absorbance of 1.0 at 260 nm in 15 min when yeast RNA is hydrolyzed at 37°C and pH 7.0. The

specific activity is >2.5×106 unit/mg.

Purity: > 95 % as determined by SDS-PAGE

Endotoxin:  $< 1.0 \text{ EU/}\mu\text{g}$  of the protein as determined by the LAL method.

Lyophilized from a solution filtered through a 0.22 µm filter, containing PBS, pH 7.4. Typically,

Formulation: a mixture containing 5% to 8% trehalose, mannitol, and 0.01% Tween 80 is incorporated as a

protective agent before lyophilization.

## **Preparation and Storage**

#### Reconstitution:

A Certificate of Analysis (CoA) containing reconstitution instructions is included with the products. Please refer to the CoA for detailed information.

#### Stability & Storage:

It is recommended to store recombinant proteins at -20°C to -80°C for future use. Lyophilized powders can be stably stored for over 12 months, while liquid products can be stored for 6-12 months at -80°C. For reconstituted protein solutions, the solution can be stored at -20°C to -80°C for at least 3 months. Please avoid multiple freeze-thaw cycles and store products in aliquots.

# Shipping:

In general, Lyophilized powders are shipping with blue ice.

# **Protein Background**

RNase A, also known as ribonuclease A and RNASE1, belongs to ribonuclease A superfamily. It is a pancreatic-type of secretory ribonuclease. RNase A is a basic protein and its many positive charges are consistent with its binding to RNA (a poly-anion). More generally, RNase A is unusually polar or, rather, unusually lacking in hydrophobic groups, especially aliphatic ones. As an endonuclease, RNase A cleaves internal phosphodiester RNA bonds on

Page 1 of 2 www.targetmol.com

the 3'-side of pyrimidine bases. It prefers poly(C) as a substrate and hydrolyzes 2',3'-cyclic nucleotides, with a pH optimum near 8.0. RNase A is monomeric and more commonly acts to degrade ds-RNA over ss-RNA. Alternative splicing occurs at this locus and four transcript variants encoding the same protein have been identified.

### Reference

Tubert P,et al. (2011) Interactions crucial for three-dimensional domain swapping in the HP-RNase variant PM8. Biophys J. 101(2):459-67.

Vinayagam A,et al. (2011) A directed protein interaction network for investigating intracellular signal transduction. Sci Signal. 4(189):rs8.

Fischer S,et al. (2011) Expression and localisation of vascular ribonucleases in endothelial cells. Thromb Haemost. 105(2):345-55.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481

Page 2 of 2 www.targetmol.com