

#### **TECHNICAL DATA SHEET**

Equivalent Performance, Exceptional Value

# **Cell Stimulation Cocktail (500X)**

Catalog Number: TNB-4975

## PRODUCT INFORMATION

**Contents:** Cell Stimulation Cocktail (500X)

Use By: 6 months from date of receipt


Storage Conditions: -20°C, protected from light and

tightly sealed to prevent evaporation

Formulation: 500X concentrate of Phorbol 12-

> Myristate 13-Acetate, Ionomycin, Brefeldin A, and Monensin in

Ethanol



C578l/6 splenocytes were stimulated with 1X Cell Stimulation Cocktail (500X) (TNB -4975-UL100) (right panel) or unstimulated (left panel) and then stained with APC Anti-Mouse CD3e (20-0031), followed by intracellular staining with PE Anti-Mouse IL-2 (50-2031).

#### DESCRIPTION

Cell Stimulation Cocktail (500X) is a ready to use cocktail that can be used for polyclonal cell activation. It is composed of the phorbol ester, PMA (Phorbol 12-Myristate 13-Acetate), a calcium ionophore (Ionomycin), and the protein transport inhibitors Brefeldin A and Monensin. This Cell Stimulation Cocktail (500X) can be used to elicit cytokine production and retention in a variety of cell types allowing for detection of intracellular protein by flow cytometry.

#### **PREPARATION & STORAGE**

Cell Stimulation Cocktail is supplied as a 500X solution in Ethanol. Product should be stored at -20°C, protected from light and tightly sealed to prevent evaporation.

### **APPLICATION NOTES**

Cell Stimulation Cocktail (500X) has been quality-tested by treatment of mouse splenocytes for 6 to 18 hours followed by intracellular staining and flow cytometric analysis of IL-2 or Ifn-gamma. Cell Stimulation Cocktail is supplied as a 500X solution and should be used at 1X final concentration for cell stimulation. It is recommended to add the cocktail directly to the cell suspension at 2 uL/mL to obtain a 1X final concentration.

Tonbo Biosciences tests all antibodies by flow cytometry. Citations are provided as a resource for additional applications that have not been validated by Tonbo Biosciences. Please choose the appropriate format for each application and consult Materials and Methods sections for additional details about the use of any product in these publications.

For Research Use Only

Not for use in diagnostic or therapeutic procedures. Not for resale. Not for distribution without written consent. Tonbo Biosciences will not be held responsible for patent infringement or other violations that may occur with the use of our products. Tonbo Biosciences, Tonbo Biosciences Logo and all other trademarks are the property of Tonbo Biotechnologies Corporation. © 2013 Tonbo Biosciences.